



# MGC Health™ D-PEN MINI Marijuana Drugs of abuse Saliva Oral Test 2 Tests/Kit

MGC Health™ D-PEN MINI Marijuana DRUG of Abuse Saliva Oral test is a rapid screening test. It's a lateral flow, one-step immunoassay for the quality detection of cannabinoids in human oral fluid at the cut-off level of 40 ng/mL

## Key Product Features

1. Advanced technology kits for precise testing for Drugs of Abuse.
2. Tamper resistant tests provide results in 5 minutes
3. Full observation during collection allows test to be performed anytime and anywhere.
4. FDA 510(k) exempt for Employment & Insurance Use



## Material Provided

1. D-Pen Mini THC Oral Fluid Test
2. Instruction for use

## Storage and stability

1. Store at 35°F - 86°F (2°C - 30°C) in the sealed pouch up to the expiration date
2. Keep away from direct sunlight, moisture and heat
3. DO NOT FREEZE
4. Preferably open the pouch only shortly before the test

## Test Procedure and Interpretation



## Ordering Information

Specification	Description
UPC	383173000252
MGC Ref	IT0000106
Retail Box Shelf dimensions	5.9 x 2.4 x 1
Shipping Dimensions	11.9 x 7.5 x 4.5
Units per Case	24
Cases per pallet	40

**Available Display for  
Drug of Abuse tests**

